

signia

Life sounds brilliant.



The Signia hearing care professional support - Overview



Hearing
Systems

SIEMENS

Audio streamers.

Products	Type	Description	Item No.
StreamLine TV 	StreamLine TV	StreamLine TV connects TVs without Bluetooth functionality to Signia Nx devices and Pure 13 BT primax hearing aids for wireless audio streaming.	10944539
StreamLine Mic 	StreamLine Mic New!	The new StreamLine Mic provides full hands-free audio streaming of phone calls and music with any Bluetooth smartphone. It also works as a remote microphone and is compatible with Signia Nx hearing aids.	10945010
easyTek and easyTek App 	easyTek	easyTek audio streamer and remote control incl. neckloop antenna. Transmitters and VoiceLink need to be purchased separately. Incl. silicone sleeve for protection against moisture. Size: 49x12 mm	10931133
	easyTek App	The easyTek App allows Android and iOS users to control their easyTek via their smartphone.	Free download
TV Transmitter 	TV Transmitter	The TV Transmitter connects TVs without Bluetooth functionality to the easyTek for wireless audio streaming.	10937963
VoiceLink 	VoiceLink set	Transmitter, lapel microphone, power supply set, audio cable set, adapter set	10673506
	Microphone service set	Lapel microphone, accessory set	10673511
CROS Pure 	CROS Pure	For people with unaidable hearing loss on one side. CROS Pure receives sound on that side and transmits it to an hearing aid on the other side. Compatible with almost all primax hearing aids.	10934252 (Beige) 10941192 (Brown) 10945867 (Dark Brown) 10941196 (Deep Red) 10941193 (Granite) 10941194 (Grey) 10941195 (Silver)

Remote controls and eCharger.

Products	Type	Description	Item No.
<p>myControl App</p> 	myControl App	Compatible with Android and iOS, the myControl App allows wearers of Signia Nx devices and Pure 13 BT primax to enjoy the best hearing experience via Bluetooth. It also offers independent volume control, the ability to select which notification tones they hear, sound exposure and wear time tracking.	Free download
<p>touchControl App</p> 	touchControl App	The touchControl App allows Android and iOS users to control their hearing instrument programs, volume, bass and treble via their smartphone.	Free download
<p>easyPocket</p> 	easyPocket	User-friendly, ergonomic remote control with a modern design and built-in display. Size: 92.5x42.5x17.5 mm	10675287
<p>miniPocket</p> 	miniPocket	Small enough to carry on a keyring, this discreet, convenient remote control avoids the need for a large remote control or a smartphone.	10939860
<p>eCharger</p> 	<p>eCharger with electronic drying function</p> <p>Inlay 13-K</p> <p>Inlay 13-M</p> <p>Inlay 13-N</p> <p>Inlay 312-L</p> <p>Rechargeable battery 312</p> <p>Rechargeable battery 13</p> <p>Car adapter</p> <p>KuanTech power supply</p>	<p>eCharger without inlay. Rechargeable battery not included.</p> <p>Inlay for Motion PX binax/micon</p> <p>Inlay for Carat binax</p> <p>Inlay for Motion SX primax/binax</p> <p>Inlay for Pure primax/binax/micon</p> <p>For size 312 rechargeable hearing aids</p> <p>For size 13 rechargeable hearing aids</p> <p>For eCharger, miniTek and Tek</p> <p>With 4 power supply adapters</p>	<p>10600230</p> <p>10823023</p> <p>10824799</p> <p>10827814</p> <p>10820446</p> <p>10423863</p> <p>10423862</p> <p>10601544</p> <p>10666693</p>

miniReceiver 2.0

More robust, more flexible.

Compatible hearing aids

- Ace primax	- Ace binax
- Pure 13 BT primax	- Pure binax
- Pure primax	- Carat binax
- Carat primax	- Orion 2 RIC
- Cellion primax	- Orion RIC

The miniReceiver 2.0 is a smart RIC receiver that uses SIM technology. It securely and reliably stores and transmits data between the receiver and the hearing aid. The SIM card can also detect an incorrect receiver and automatically reduce gain to protect the wearer. And the design enables more ergonomic placement of the receiver for better microphone placement.

miniReceiver 2.0 - ready to perform each and every day.

Robust connector
Protects against dust, dirt, and moisture.

Color indicator
Blue left side and red right side immediately identifiable.

Reliable click interface
Compatible with the premium Click Sleeves for a comfortable, immediate fit, the interface also offers a proven connection to Click Domes™ and Click Molds™.

Effective 2-step wax protection
In combination with Click Domes and Click Molds.


SIM Card Smart technology

- Auto identification safety feature
- Reliable and secure logical connection
- Enables better sound quality

- 1 Front microphone moved to top
- 2 SIM card
- 3 Plug and jack moved down

HF4 Pro filter
The high-tech filter gauze reliably protects the hearing instrument by repelling dirt, cerumen, oil and water - for an extra-long usability.

miniReceiver 2.0 - exchanging is as easy as 1, 2, 3.



Open the locking bar.


Remove/Insert the miniReceiver 2.0.

Close the locking bar.

The miniReceiver 2.0 is easily exchangeable. To replace the miniReceiver 2.0, use the exchange tool to open and close the locking bar for reliable securing and no risk of dislodgement. To see an animation of how to exchange the miniReceiver 2.0 visit the website for professionals.

miniReceiver fitting.

Available in four power forms and in different lengths, our broad range of miniReceivers caters to nearly all your customers' listening requirements and wearing preferences.



Products

miniReceiver 2.0 Fitting set

Content

miniReceiver S/M/P each with size 2 left and right; Click Dome packs each with 8 mm open, 8 mm closed, 8 mm double, semi-open; Click Mold service set, filter HF4 left and right, set of miniReceiver wax guards, measuring card for receiver length, measuring card for Click Mold size

Item No.

10670779

Type S



Size

Units/box

miniReceiver 2.0

miniReceiver

0L	1	10669907	10425988
0R	1	10669908	10425989
1L	1	10669909	10425990
1R	1	10669910	10425991
2L	1	10669911	10425992
2R	1	10669912	10425993
3L	1	10669913	10425994
3R	1	10669914	10425995

Type M



0L	1	10824820	10671603
0R	1	10824821	10671604
1L	1	10824822	10671605
1R	1	10824823	10671606
2L	1	10824824	10671607
2R	1	10824825	10671608
3L	1	10824826	10671616
3R	1	10824827	10671619
4L	1	10933054	-
4R	1	10933053	-

Type P



0L	1	10669923	10426004
0R	1	10669924	10426005
1L	1	10669925	10426006
1R	1	10669926	10426007
2L	1	10669927	10426008
2R	1	10669928	10426009
3L	1	10669929	10426010
3R	1	10669930	10426011








Type HP



Please order the miniReceiver size HP with a custom CIC shell in the same way as a custom hearing aid using the order form.


* Only for miniReceiver 2.0 M.

Click Domes and Click Molds.

Products	Description	Size	Units/box	Item No.
Click Dome Types 	<p>The open dome is for patients who require maximum comfort and no occlusion.</p>	4 mm	6	10426024
		8 mm	6	10426022
		10 mm	6	10426023
	<p>The semi-open dome combines the benefits of an open and a closed dome for patients who desire the comfort of an open fitting but need more power.</p>	8-12 mm	6	10669971
	<p>The closed dome is for patients who require more loudness than the open dome can offer.</p>	6 mm	6	10426025
		8 mm	6	10426020
		10 mm	6	10426021
	<p>The double dome is for patients who require broadband amplification and more retention. Can also be used with a miniReceiver size P.</p>	8/10 mm	6	10426026
		10/12 mm	6	10426027
Click Mold Types  <p>Short mold Open mold, short</p> <p>Long mold Open mold, long</p>	<p>Click Molds are an alternative to fitting with Click Domes. These molds are made by a certified partner. If you have any questions, please do not hesitate to contact us.</p>			
Click Mold User Set	<p>Short Click Mold disassemble tool, waxpick and magnet, receiver mold tube set small</p>			10297177
	<p>Long Click Mold disassemble tool, receiver long mold disassemble set</p>			10420833
Click Mold Starter Kit 	<p>Click Mold Starter Kit Package contents: Click Mold adapters (red and blue), Click Mold service set, Click Mold disassembly tool, Click Mold assembly tool, long mold drilling tool, molding cutters, removal strings, miniReceiver S/M/P dummies, Click Mold measuring card, HF4 filters left and right, miniReceiver wax guard set, work instructions on CD</p>			10426787
	<p>Click Mold Starter Kit for miniReceiver 2.0</p>			10670820
	<p>Click Mold measuring gauge</p>			10598784
	<p>Click Mold measuring gauge for miniReceiver 2.0</p>			10670821

Products

Accessories




Description

Item No.

1. HF4 Pro filters right (red) (10 pcs) and left (blue) (10 pcs)	10827324
2. HF4 filters left (blue) (20 pcs)	10600985
3. HF4 filters right (red) (20 pcs)	10600984
miniReceiver wax guard (10 pcs)	10597971
1. miniReceiver 2.0 removal tool	10674370
2. Click Dome removal tool	10605850
Click Mold service set	10597930
Template for receiver length	10174571
miniReceiver 2.0 measuring gauge For receiver length & impression measurement for P/HP receivers	10673959
miniReceiver HP measuring gauge For receiver length & impression measurement for HP receivers	10670257
Nitro micon measuring gauge The measuring gauge supports the fitting process by letting you easily find out about the right size for the wearer's custom hearing aids.	10823745

New! Click Sleeves 2.0.

One click and it fits. The exchangeable Click Sleeves 2.0 can easily be clicked into place. Made from super-soft silicone and smoky-colored material, they are available in four different sizes that fit in almost every ear canal for a simple, quick and swift fitting. In addition, turn protection and side indication provide an extraordinary wearing comfort.



Click Sleeves 2.0 are available for all Signia RIC and Click ITE hearing aids.

Products

Content

Item No.

Type

Size

vented



XS (6 pcs)	10951144
S (6 pcs)	10951041
M (6 pcs)	10951043
L (6 pcs)	10951045

closed



XS (6 pcs)	10951145
S (6 pcs)	10951042
M (6 pcs)	10951044
L (6 pcs)	10951046


ThinTubes.

Slim shape, secure fit, superb performance.

A secure fit and a simple exchange process in just a few easy steps: ThinTubes have a threaded screw connector that is easy to attach and remove from the housings of Motion SX/SA and P primax and Motion SX/SA binax models. In combination with the ThinTube Adapter, it is also compatible with Motion S primax and a range of other hearing aid models.

Connector

Firmly attached to or integrated into the housing. The interface with a screw thread ensures a secure fit between the hearing aid and the ThinTube or the earhook.



Thanks to the interface with a screw thread, the ThinTube can be quickly and easily screwed on or off.

Attaching the ThinTube.



1 Move the ThinTube close to the interface on the hearing aid. Make sure the ThinTube is straight on the screw thread.

2 Rotate the ThinTube clockwise until it sits fully tightened in the screw thread.


3 Stop when the ThinTube sits firmly in the interface.

4 Please check if the ThinTube is aligned to the hearing aid body to ensure a perfect fit and optimum sound transmission.



Designed for a high cosmetic appeal:

A slimmer design and a more transparent material make the ThinTube more discreet and attractive than comparable solutions.



Motion S

To make Motion S binax and primax models compatible with a ThinTube, their standard interface can be easily exchanged by screwing a ThinTube Adapter on.

Product

Content

Item No.


ThinTube Adapter



Screw on adapter for all products with LifeTube to make them compatible for usage with all ThinTubes. Compatible hearing aids: Motion S primax, Orion 2, Sirion 2, Intuis 2, Orion, Intuis

10932441

Removing the ThinTube.



To remove the ThinTube, rotate it counterclockwise until it is fully detached from the housing.

Acoustic coupling for ear molds.



2 in 1 earhook
Compatible with Motion SX/SA primax and binax.

Item No. 10935294



Earhook
Compatible with all ThinTube ready devices.

Earhook Item No. 10828811
Earhook small Item No. 10933628

LifeTube, S-LifeTube, ThinTube.

Products

LifeTube

Item No.

Compatible BTE hearing aids: Motion S and P primax,
Orion 2 S/M/P, Sirion 2 S/M/P

Fitting Set

Content



LifeTubes in all lengths, LifeTips in all sizes,
LifeTube accessory box incl. measuring gauge,
LifeTube shaper, cleaning wires

10296968

LifeTube, S-LifeTube and ThinTube Refill Boxes



LifeTube

S-LifeTube

ThinTube

Type

Units/box

1L	10	10296944
1R	10	10296945
2L	10	10296946
2R	10	10296947
3L	10	10296948
3R	10	10296949
4L	10	10296950
4R	10	10296951
5L	10	10296952
5R	10	10296953
6L	10	10296954
6R	10	10296955


LifeTips.

Products

Content

Item No.

LifeTip Refill Boxes



Open tips

Closed tip

Double tips

Semi-open tip

Type

Units/box

open 4 mm	10	10296987
open 6 mm	10	10296988
open 8 mm	10	10296989
open 10 mm	10	10296990
closed 8 mm	10	10296986
double 8/10 mm	10	10296984
double 10/12 mm	10	10296985
semi-open 8-12 mm	10	10673779

Additional Accessories

Length

Units/box

LifeTube shaper



1	5	10054903
2	5	10054904
3	5	10054905
4	5	10148704
-	30	10054970
-	1	10054711

Life cleaning wires

Measuring gauge

S-LifeTube	Item No.	ThinTube	Item No.
Compatible BTE hearing aids: Motion PX binax, Motion M/P/PX micon		Compatible BTE hearing aids: Motion SX, SA and P primax, Motion SX/SA binax	
Content		Content	
S-LifeTubes in all lengths, LifeTips in all sizes, LifeTube accessory box incl. measuring gauge, LifeTube shaper, cleaning wires	10296992	ThinTubes in all lengths, measuring gauge, cleaning wires	10828299


Type	Units/box	
1L	10	10297010
1R	10	10297011
2L	10	10297012
2R	10	10297013
3L	10	10297014
3R	10	10297015
4L	10	10297016
4R	10	10297017

Type	Units/box	
1L	5	10828291
1R	5	10828292
2L	5	10828293
2R	5	10828294
3L	5	10828295
3R	5	10828296
4L	5	10828297
4R	5	10828298
5 side neutral	5	10828971

Lotus Coupling.

Products	Content	Item No.
Type		
Comfort coupling set	For Prompt P and SP, Run, Fun	10937494
Simple coupling set	For Fast	10937958

Earmolds.

Products	Content
Earmolds for classic fitting	
	We offer a comprehensive range of acrylic earmolds for all BTE models. Selected styles are also available as soft earmolds: Skeleton, Canal, Full shell, Half shell.
<p>Skeleton Half helix Full helix Tube holder</p>	

iScan II.

Smaller. Faster. More accurate.




Invaluable for ensuring a perfect fit on the spot, iScan II lets you scan and view 3D images of impressions for ITE instruments or ear molds. Smaller, approximately 50% faster, and even more accurate, iScan II saves you time and money. Its increased reliability means that accuracy checks are a thing of the past. Making adjustments and, if necessary, fine tuning the virtual 3D model is a snap, reducing the number of impressions needed per customer. All it takes is a few clicks and the hearing instrument is ordered online, reducing processing time by 48 hours. Plus, you can track the entire process online from ordering and manufacturing to delivery.



iScan II benefits at a glance:

- Compact design
- Larger operating space
- Accurate to within 1/50th of a millimeter
- Compatible with Windows 7 to 10
- Compatible with Noah 3 and 4
- Standalone capability
- Quick and safe ordering via e-mail
- Open data format

Three steps: From 3D to despatch.



1. Scan

Secure the ear impression to the scanning arm in the iScan II and start the program. Thanks to the high resolution scanner technology, the scan results are accurate to within 1/50th of a millimeter.

2. Check

In about 60 to 70 seconds, the precisely reconstructed impression will appear on the screen as a 3D scan. Now you can check it from all angles and, if necessary, use the software to add comments for making adjustments.

3. Order

A few clicks are all it takes and your order is on its way to us by e-mail. Your orders are processed simply, efficiently and, from a time perspective, completely flexibly. Plus, you can track the complete process online.



Type

iScan II

Product details	Ear impression scanner	
Dimensions	Approx. L 341xW 320xH 125 mm	
Weight	Approx. 7 kg	
Technology	Stripe-light scanner with durable high-performance LED	
Number of impressions per scan process	1	
Supply voltage	100 - 240 VAC; 50/60 Hz	
Power consumption	Max. 30 W	
Protection category	IP11	
Accuracy	< 20 µm	
Duration of scan	Approx. 60 - 70 sec. when following recommended computer performance	
Output data format	ASCII or STL	
Ports	USB 2.0; USB 3.0	
Scope of delivery	Scanner, power cable, USB cable, CD with iScan II software, calibration data and user manual, 1 impression holder, printed version of Quick Installation Guide	
System requirements	Minimum	Recommended
System	Windows 7	Windows 10
Ram	2 GB RAM	4 GB RAM
USB	USB Port 2.0	USB Port 3.0
Graphic card	Graphic card minimum 256 MB RAM	Graphic card 512 MB RAM

Products

Description

Item No.

iScan II



iScan II

- Larger operating space
- Optimized scanning process: Up to twice as fast as the previous version, approx. 60 - 70 sec.
- Encrypted and open scanning mode
- Works with all impression colors; transparent material needs to be colored with a spray

10936441

Software

- Full integration options: Compatible with Windows 7 to 10, Noah 3 and 4, or standalone capability

10936443



The soundDemo Toolkit:
This unique set of valuable tools makes your fitting workflow smoother, more seamless and even more successful.

soundDemo Toolkit.

Demonstrating primax benefits to your customers – from pre-sales to home trials and in your office.

The soundDemo Toolkit contains a range of unique tools that support your customer acquisition activities and help smooth your workflow throughout the entire fitting process. soundDemo Videos, soundDemo Suite, and soundDemo Hearing Aids fit seamlessly in your existing workflow, saving you time and effort during the awareness, orientation, fitting and purchase phases.

Plus, they help you turn the process into a premium experience for your customers that's both informative and entertaining.

Three phases of the customer experience:

AWARENESS

soundDemo Videos

are the perfect tool for helping sufferers understand hearing loss and the difference that primax hearing aids can make.

ORIENTATION

soundDemo Suite

allows customers to experience different auditory situations through hearing aids at first hand, bringing the benefits of the primax platform to life.

FITTING & PURCHASE

soundDemo Hearing Aids

can be reconfigured as 7px, 5px or 3px instruments so that only one pair is needed for the entire trial period, eliminating the need for multiple devices.

Once the trial is over, all of the customer's individual settings can be transferred from the soundDemo Hearing Aids to the customer's own, personalized instruments quickly and simply via Connex Eight.

Reduced effort and smoother workflow:


The **soundDemo Toolkit** saves time and reduces costs throughout the fitting process.

soundDemo Videos.

Using sound recordings taken directly from hearing aids, soundDemo Videos show real life situations and demonstrate - visually and audibly - the benefits of primax features and programs in a way that's easy to understand.


Ideal for use before, as well as during consultations, they can be viewed and experienced online to create greater awareness of the benefits of primax hearing aids and help you reach new customers.

Find all soundDemo Videos at [youtube.com/c/signia-hearing](https://www.youtube.com/c/signia-hearing)



Objectively recorded:


A dummy fitted with Pure primax hearing aids and a headset camera is used to record real-life situations that show the performance of primax hearing aids.



soundDemo Suite.

The ideal tool during the fitting process, soundDemo Suite gives an unbiased demonstration of the benefits of primax hearing aids in the most convincing way possible. Using either primax or soundDemo Hearing

Aids after First Fit, speech, music and background sounds are played via a DVD and 5.1 surround sound system for a truly immersive experience.



Every hearing situation in one place:

Various hearing situations clearly demonstrate the features of primax hearing aids, e.g. while driving, in a restaurant or listening to music.



soundDemo Hearing Aids.


Perfect for live demonstrations, soundDemo Hearing Aids are designed for fitting and home trials, allowing your customers to experience their full potential. Based on primax hearing aids, they are available as RIC, BTE and Click CIC models and function the same way as normal primax hearing instruments.

The timer function in soundDemo Hearing Aids lets you define the length of the home trial to suit your and your customer's needs - from one to six weeks in weekly steps. At the end of the trial period, a repeating signal is heard reminding the customer to return the trial hearing aid.




One pair of instruments, a variety of options.

All soundDemo Hearing Aids can be configured as 7px, 5px or 3px instruments, reducing stock complexity. They are ideal for letting customers try out accessories for themselves and can also be used to provide wearers with temporary replacement instruments if a hearing aid needs to be repaired.



soundDemo Hearing Aids are recognizable by the "sDemo" identifier on the product and in Connex Eight.

soundDemo Hearing Aids:
One model, all performance levels.




Connex Eight: Programming soundDemo Hearing Aids.

The perfect partner for configuring soundDemo Hearing Aids, Connex Eight lets you quickly and conveniently adjust and store individual settings throughout the fitting, even when changing performance levels. At the end of the process, you can effortlessly transfer the settings to the customer's own hearing aids for a streamlined workflow.



When a soundDemo instrument is connected, an icon appears that takes you to the special soundDemo configuration page.

Performance level selection



Timer adjustment for setting trial duration

Products

Description

Item No.

Ace sDemo

Light and ultra-discreet RIC instrument in a compact form

10934292

Pure 13 BT sDemo

Elegant RIC instrument providing connectivity beyond streaming

10941118

Pure sDemo

Fully featured binaural RIC instrument for discerning wearers

10936013

Carat sDemo

Elegant and powerful RIC instrument with a high level of amplification

10944638

Motion SX sDemo

Sleek and robust BTE instrument with rechargeable battery

10936219

Motion SA sDemo

Convenient and comfortable BTE instrument with direct audio input option

10934923

Motion P sDemo

Powerful and small BTE instrument for severe to profound hearing losses

10936891

Silk sDemo left

Swift fit, super-discreet blue Click CIC for the left side



10938585

Silk sDemo right

Swift fit, super-discreet red Click CIC for the right side

10938584

Unity 3.

Products	Description	Item No.
Unity 3 Fitting Unit System 	Unity™ 3 Fitting unit package incl. Unity 3 fitting unit, DC adapter, USB cord, client response, Unity quick installation guide, Unity software USB stick	10823931
Unity 3 HIT package 	Unity 3 HIT package incl. coupler microphone, reference microphone, battery pills (S size), BTE tubes (3 pcs), ITE paddy, BTE coupler attachment, ITE coupler attachment, 2cc coupler attachment, body coupler attachment, USB cord, HI-PRO extension cable, Unity quick installation guide, Unity software USB stick	10823932

Products	Item No.	Products	Item No.
Transducers		Accessories	
Real ear probe	10823933	USB isolation cord	10823939
EAR 3A	10823934	Free field speaker (one)	10824086
DD45	10823935	Headset with mic (talk over)	10823940
Inserts earphone	10823936	Microphone (talk back)	10823941
HDA200	10823937	Client response	10823942
Bone conductor - B71	10823938	Ear hanger (4 pcs)	10823943
Consumables		DC adaptor	10823944
HIT ITE paddy	08786084	USB cord	10824087
Insert tube with nipples	10823956	Sound chamber cables	10823946
Probe tube guide (8 pcs)	10823957	HIT ITE coupler attachment	10823948
Probe tubes (36 pcs)	10823958	HIT BTE coupler attachment	10823949
BTE tube (3 pcs)	10823959	HIT 2cc coupler attachment	10823947
ITE foam	10824088	HIT body couple attachment	10823950
BTE foam	10824089	HIT coupler microphone	10823951
Battery pills		HIT reference microphone (calibrated)	10823952
Battery pills (5 sizes)	10823962	HIT HI-PRO extension cable	10824187
Battery pill 5A (1 pc)	10823967		
Battery pill 10A (1 pc)	10823968		
Battery pill 312A (1 pc)	10823969		
Battery pill 13A (1 pc)	10823970		
Battery pill 675A (1 pc)	10823971		

Workflow and solutions.

Products	Description	Item No.
Wireless Programming  <p>ConnexxAir Noahlink Wireless</p>	ConnexxAir Convenient wireless programming system for hearing aids equipped with e2e wireless capabilities that are programmable in Connexx 8. Not compatible with Pure 13 BT primax. Includes a neck loop, programmer, cradle and Bluetooth dongle	10944827
Programming  <p>HI-PRO 2 USB Noahlink</p>	HI-PRO 2 USB PC Interface Hearing aid interface with HI-PRO cable, using the serial USB port for communication between a PC and the HI-PRO USB unit HI-PRO 2 cable right HI-PRO 2 cable left Noahlink incl. Bluetooth adapter Noahlink , wireless Bluetooth communication device to program hearing aids with ease, both in office and on location Noahlink cable	10297179 02944986 02944994 07784890 L: 07784908 R: 07784916
Programming Adapter 	Programming adapter, size 10A (1 pc) Programming adapter, size 312 (1 pc) Programming adapter, size 13 (1 pc) Programming adapter, Flex for custom (1 pc)	10824471 10824469 10824470 03726994
Programming Shoe	Programming shoe for Sirion 2 SP, Intuis 2 SP and Prompt SP Programming shoe for Prompt P	10940144 10933080
soundDemo Suite DVD 	The integrated system for demonstrating hearing aid features. All you need is a <ul style="list-style-type: none"> - 5.1 surround loudspeaker set up - pair of primax instruments or soundDemo Hearing Aids. Available in English, German, Japanese, Chinese, Korean, Dutch, Portuguese.	PAL version: 02991-99E1-X1 NTSC version: 03020-99E1-X1
Counseling Suite 	Counseling Suite The comprehensive advisory tool with practical functions. <ul style="list-style-type: none"> - Helpful support during the consultation process with customers - before, during and after the fitting - Assistance in demonstrating and familiarizing customers with the topics of "hearing" and "hearing instruments," e.g. the hearing loss simulator, selection of the right performance level, maintenance tips - With content on tinnitus, binaural hearing and improved upselling support - Intuitive design for easy use 	Available for download at signiapro.to/counselingsuite at no cost.
	Counseling Suite App	Available on the Google Play Store or Apple App Store at no cost.

Legal Manufacturer

For Siemens products:
Sivantos GmbH
Henri-Dunant-Strasse 100
91058 Erlangen
Germany

Sivantos GmbH is a Trademark
Licensee of Siemens AG.

For Signia products:
Signia GmbH
Henri-Dunant-Strasse 100
91058 Erlangen
Germany

Signia GmbH is a Trademark
Licensee of Siemens AG.

© 11.2017 Signia GmbH

Find us on:
[signia-hearing.com](https://www.signia-hearing.com)



facebook.com/signiahearing



youtube.com/signiahearing